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
AllCheck COVID19 Ag

※ Please read the instructions for use carefully.

[INTENDED USE]

AllCheck COVID19 Ag is a rapid immunochromatographic test for the qualitative detection of specific antigens to the SARS-CoV-2 in human nasopharynx.

[SUMMARY]

Coronavirus disease 2019 (COVID-19) is an infectious disease caused by SARS-CoV-2, a new strain of coronavirus that has not been previously identified in humans. The disease is primarily spread between people via respiratory droplets from infected individuals. Time from exposure to onset of symptoms is generally between 2 and 14 days. The disease may initially present with few or no symptoms, or may develop into fever, coughing, shortness of breath, pain in the muscles and tiredness.

[TEST PRINCIPLE]

The “**AllCheck COVID19 Ag**” is an in vitro diagnostic tests that measures the presence of COVID-19 Ag(SARS-CoV-2 nucleocapsid antigen) in a person’s nasopharyngeal swab using the immunochromatographic assay and is used as a test for diagnosing coronavirus infection.

The test has two pre-coated lines, a Control line (C) and a Test line (T), on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimens. When the specimen is added to sample pad, it moves through the conjugate pad and mobilizes gold conjugates: COVID-19 Ag (Nucleocapsid protein) specific antibody conjugate that is coated on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with anti-COVID-19 Ag that is coated on the test region. If COVID-19 antigen is present, the result is the formation of a colored band in the test region. If there is no COVID-19 antigen in the sample, the area will remain colorless. The sample continues to move to the control area and forms a red or purple color, indicating the test is working and the result is valid. This test is intended for professional use as an aid in the diagnosis of COVID-19 antigen.

[REAGENTS AND MATERIAS]

- COVID19 Ag Test devices
 - Extraction buffer tube
 - Nozzle cap
 - Sterile swab
 - Instructions for use
- Materials not provided but required: Gloves, Timer.

[STORAGE AND STABILITY]

AllCheck COVID19 Ag can be stored for up to 18 months when stored at between 1~30℃(34~86°F). Keep the test cassette in sealed pouch until use. Do not open the pouch until you are ready to perform a test. Do not use the test after the indicated expiration date.

[WARNINGS AND PRECAUTIONS]

- For In Vitro Diagnostic Use only. Do not re-use test device.
- The product should only be used by trained clinical professionals.
- Do not use the device past its expiration date.
- To obtain accurate results, you must follow the package insert instruction.
- Store the devices at 1~30℃(34~86°F). Do not freeze.
- If you store the product at refrigerator, leave it at room temperature before testing.
- Do not use the test kit if its packing is damaged or the seal is broken.
- Use the devices immediately after opening the aluminum pouch.
- Do not touch the membrane zone.
- Do not mix the components from different lot to lot.
- Do not eat or smoke while handling specimens.
- Wear protective gloves and mask while handling specimens. Wash hands thoroughly afterwards.
- Avoid splashing or aerosol formation.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container, as if they were infectious waste.
- All waste except liquid should be autoclaved at 121℃ for more than 1 hour or decontaminated and disposed in accordance with local and international regulations.
- Please calibrate the pipette when you transfer the samples.
- SARS-CoV-2 virus infection cannot be diagnosed by the results of this device alone.
 - False positive test results may occur due to cross-reactivity. Even if a positive result is obtained from the device, the final confirmation of infection should be made by an expert by checking clinical symptoms, and using other test methods.
 - When the antigen concentration is low, a negative result can occur. The possibility of infection cannot be ruled out even if the device indicates a negative result.
- This device only detects the presence of SARS-CoV-2 nucleocapsid protein antigen, and there is no correlation between the strength of the test line.

[SPECIMEN COLLECTION AND PREPARATION]

- To collect a nasopharyngeal swab specimen, insert a sterile swab into the nostril of the patient, reaching the surface of the posterior nasopharynx.
- Using gentle rotation, push the swab until resistance is met at the level of the turbinate.
- Rotate the swab a few times against the nasopharyngeal wall.
- Remove the swab from the nostril carefully.
- Specimen should be tested as soon as possible after collection.
- Use the collected specimen and extraction buffer immediately. Be careful of contamination.
- Specimens may be stored at room temperature for up to 1 hours or at 2-8℃(36-46°F) for up to 4 hours prior to testing.
- Transport fresh specimens to the laboratory as rapidly as possible in a suitable liquid transport system.
- For nasopharyngeal swabs in transport media, use the 1-3 ml of transport media.

[LIMITATIONS]

- The test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
- False negative results may occur if a specimen is improperly collected, transported, or handled. Specimens should be test as quickly as possible after specimen collection.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other non-SARS viral or bacterial infections.

A negative result does not preclude the possibility of infection with COVID-19. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.

[PERFORMANCE CHARATERISTICS]

- Clinical performance evaluation

In order to test the detection sensitivity and specificity of this test, 100 nasopharyngeal swab specimen samples (40 Positives and 60 Negatives) were collected from individuals who tested with a RT-PCR method for SARS-CoV-2 infection. 100 nasopharyngeal swab specimen samples were tested with **AllCheck COVID19 Ag** and the results were summarized in the table below:

Clinical performance		RT-PCR method		Total
Result		Positive	Negative	
AllCheck COVID19 Ag	Positive	39	0	39
	Negative	1	60	61
Total		40	60	100

Sensitivity: 97.5%(95% CI: 97.12~99.56%), Specificity: 100%(95% CI: 93.98~100%)

- Limit of Detection (LoD)

The Limit of Detection (LoD) of the **AllCheck COVID19 Ag** was determined using limiting dilutions of heat-inactivated SARS-CoV-2, isolate USA-WA1/2020. The LoD is 2.5 X 10² TCID₅₀/ml.

- Cross Reactivity

Cross-reactivity of the **AllCheck COVID19 Ag** was evaluated by testing a panel of high prevalence respiratory pathogens that could potentially cross-react with the **AllCheck COVID19 Ag**. Each organism and virus was tested in triplicate. The final concentration of each organism is documented in the following table.

Cross reacting substance	Concentration	Cross-Reactivity
Influenza, A/Beijing/262/95 (H1N1)	4000ngHA/ml	No
Influenza, A/Brisbane/02/2018 (H1N1)	4800ngHA/ml	No
Influenza, A/Brisbane/59/2007 (H1N1)	8300ngHA/ml	No
Influenza, A/California/7/09 (H1N1)	4500ngHA/ml	No
Influenza, A/Canada/6294/09 (H1N1)	1x10 ^{4.23} TCID ₅₀ /ml	No
Influenza, A/Denver/1/57 (H1N1)	5x10 ^{7.25} CEID ₅₀ /ml	No
Influenza, A/Michigan/45/2015 (H1N1)	4800ngHA/ml	No
Influenza, A/New Caledonia/20/99 (H1N1)	6500ngHA/ml	No
Influenza, A/New Jersey/8/76 (H1N1)	5x10 ^{6.25} CEID ₅₀ /ml	No
Influenza, A/NWS/33 (H1N1)	5x10 ^{6.17} CEID ₅₀ /ml	No
Influenza, A/PR/8/34 (H1N1)	1x10 ^{4.62} TCID ₅₀ /ml	No
Influenza, A/Singapore/63/04 (H1N1)	1x10 ^{4.39} TCID ₅₀ /ml	No
Influenza, A/Solomon Islands/3/2006 (H1N1)	5700ngHA/ml	No
Influenza, A/Taiwan/42/06 (H1N1)	1x10 ^{4.62} TCID ₅₀ /ml	No
Influenza, A/NY/01/09 (H1N1)	1x10 ^{6.14} TCID ₅₀ /ml	No
Influenza, A/Brisbane/10/07 (H3N2)	1x10 ^{4.07} TCID ₅₀ /ml	No
Influenza, A/Kansas/14/2017 (H3N2)	6700ngHA/ml	No
Influenza, A/Hiroshima/52/2005 (H3N2)	3200ngHA/ml	No
Influenza, A/Hong Kong/4801/2014 (H3N2)	6000ngHA/ml	No
Influenza, A/Hong Kong/8/68 (H3N2)	5x10 ^{4.5} TCID ₅₀ /ml	No
Influenza, A/New York/55/2004 (H3N2)	5500ngHA/ml	No
Influenza, A/Port Chalmers/1/73 (H3N2)	2.5x10 ⁸ CEID ₅₀ /ml	No
Influenza, A/Shangdong/9/93 (H3N2)	5000ngHA/ml	No
Influenza, A/Singapore/INFIMH-16-0019/2016 (H3N2)	5400ngHA/ml	No
Influenza, A/Texas/50/2012 (H3N2)	7400ngHA/ml	No
Influenza, A/Victoria/210/2009 (H3N2)	3600ngHA/ml	No
Influenza, A/Wisconsin/67/2005 (H3N2)	2.3x10 ⁸ CEID ₅₀ /ml	No
Influenza, A/Singapore/1/57 (H2N2)	4500ngHA/ml	No
Influenza, A/Equine/Miami/63 (H3N8)	2700ngHA/ml	No
Influenza, A/Cambodia/RO405050/2007 (H5N1)	9300ngHA/ml	No
Influenza, A/Vietnam/1194/2004 (H5N1)	6700ngHA/ml	No
Influenza, A/Anhui/1/05 (H5N1)	9900ngHA/ml	No
Influenza, A/turkey/Turkey/1/2005 (H5N1)	8000ngHA/ml	No
Influenza, A/Duck/Singapore-Q/F119-3/97 (H5N3)	3800ngHA/ml	No
Influenza, A/mallard/Netherlands/12/2000(H7N3)	11500ngHA/ml	No
Influenza, A/Anhui/1/2013 (H7N9)	9500ngHA/ml	No
Influenza, A/chick/Hong Kong/G9/1997 (H9N2)	8400ngHA/ml	No
Influenza, B/Allen/45	1.2x10 ⁴ CEID ₅₀ /ml	No
Influenza, B/Brisbane/60/08	4200ngHA/ml	No
Influenza, B/Florida/02/2006	1x10 ^{5.58} TCID ₅₀ /ml	No
Influenza, B/Guangdong/120/2000	4000ngHA/ml	No
Influenza, B/Hong Kong/5/72	8.9x10 ⁶ CEID ₅₀ /ml	No
Influenza, B/Lee/40	1x10 ^{4.39} TCID ₅₀ /ml	No
Influenza, B/Malaysia/2506/2004	5000ngHA/ml	No

AllCheck COVID19 Ag

※ Please read the instructions for use carefully.

Cross reacting substance	Concentration	Cross-Reactivity
Influenza, B/Panama/45/90	5000ngHA/ml	No
Influenza, B/Phuket/3073/2013	6000ngHA/ml	No
Influenza, B/Taiwan/2/62	5x10 ^{5.75} TCID ₅₀ /ml	No
Influenza, B/Victoria/2/87	4000ngHA/ml	No
Influenza, B/Wisconsin/1/2010	5700ngHA/ml	No
Influenza, B/Yamagata/16/88	4300ngHA/ml	No
Adeno virus type 1	1x10 ⁶ TCID ₅₀ /ml	No
Adeno virus type 2	5x10 ⁵ TCID ₅₀ /ml	No
Enterovirus	6.8x10 ⁵ TCID ₅₀ /ml	No
Measles virus	5x10 ⁶ TCID ₅₀ /ml	No
Mumps virus	5x10 ⁶ TCID ₅₀ /ml	No
Rhinovirus type 1A	4.2x10 ⁵ TCID ₅₀ /ml	No
Bordetella pertussis	1x10 ⁶ CFU/ml	No
Corynebacterium propinquum	1x10 ⁶ CFU/ml	No
Haemophilus influenzae	1x10 ⁶ CFU/ml	No
Lactobacillus fermentum	1x10 ⁶ CFU/ml	No
Legionella pneumophila	1x10 ⁶ CFU/ml	No
Neisseria meningitidis	1x10 ⁶ CFU/ml	No
Neisseria gonorrhoeae	1x10 ⁶ CFU/ml	No
Staphylococcus epidermidis	1x10 ⁶ CFU/ml	No
Streptococcus pneumoniae	1x10 ⁶ CFU/ml	No
Streptococcus pyogenes	1x10 ⁶ CFU/ml	No
Escherichia coli	1x10 ⁶ CFU/ml	No
Respiratory syncytial virus (Type A)	3x10 ⁵ TCID ₅₀ /ml	No
Respiratory syncytial virus (Type B)	3x10 ⁵ TCID ₅₀ /ml	No
Coronavirus (229E)	1x10 ^{4.5} TCID ₅₀ /ml	No
Coronavirus (OC43)	1x10 ⁵ TCID ₅₀ /ml	No
Coronavirus (NL63)	1x10 ⁴ TCID ₅₀ /ml	No
MERS-Coronavirus (Florida/USA-2_Saudi Arabia_2014)	4x10 ⁴ TCID ₅₀ /ml	No
Parainfluenza virus (Type 1)	1x10 ⁵ TCID ₅₀ /ml	No
Parainfluenza virus (Type 2)	1x10 ⁵ TCID ₅₀ /ml	No
Parainfluenza virus (Type 3)	1x10 ⁵ TCID ₅₀ /ml	No
Parainfluenza virus (Type 4)	1x10 ⁵ TCID ₅₀ /ml	No

4. Interference


The following substances, that may be naturally present in respiratory specimens or artificially introduced into the nasopharynx, were evaluated with the **AllCheck COVID19 Ag** at the concentrations listed below and were found not to affect the test performance.

Interfering Substance	Concentration	Reactivity
Whole blood	1% (v/v)	No
Mucin	1mg/mL	No
Phenylephrine	1mg/mL	No
Oxymetazoline	10ug/mL	No
Otrivin	10% (v/v)	No
Beclomethasone	0.5mg/mL	No
Dexamethasone	5mg/mL	No
Flunisolide	1ug/mL	No
Triamcinolone	5ug/mL	No
Budesonide	5ug/mL	No
Mometasone	5ug/mL	No
Fluticasone	5ug/mL	No
Menthol	1mg/mL	No
Zanamivir	1mg/mL	No
Mupirocin	5ug/mL	No
Tobramycin	1mg/mL	No

■ COLLECTION OF SPECIMEN

[Nasopharyngeal swab]


1. Insert a sterile swab into the nostril of the patient, swab over the surface of the posterior nasopharynx. Withdraw the sterile swab from the nasal cavity.
2. Insert the swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Press the nozzle cap tightly onto the tube.
4. Apply 3~4drops (about 100 μ l) of extracted specimen to the specimen well of the test device.

[Specimens in transport media]

1. Using a micropipette, collect 350ul of specimen from the collection cup or VTM.
2. Mix the specimen with an extraction buffer.
3. Press the nozzle cap tightly onto the tube.
4. Apply 3~4drops of extracted specimen to the specimen well of the test device.

■ TEST PROCEDURE

Read the result at 10minutes.
Do not interpret the test result after 15 minutes.

■ RESULT INTERPRETATION

1. Negative Result

The presence of a red line in the control region, and the complete absence of a line in the test region of the result window

2. Positive result

The presence of red lines in the test region, and in the control region of the result window.

3. Invalid result

The complete absence of a line in the control region of the result window, regardless of the test line, indicates the device has failed to perform correctly. Re-test is recommended.

■ QUALITY CONTROL

Each Test device has a built-in control. A red colored line in the detection window at the Control line can be considered an internal positive procedural control. The Control line will appear if the test procedure has been correctly performed.

CALTH Inc.

Headquarters: #321, 54, Changeop-ro, Sujeong-gu, Seongnam-si, Gyeonggi-do, 13449, Republic of Korea
TEL: 82-31-754-0320 FAX : 82-31-754-0321

Manufacturing site: #7508, 140, Beolmal-ro, Dongan-gu, Anyang-si, Gyeonggi-do, 14057, Republic of Korea
TEL: 82-31-360-0328 FAX : 82-70-8228-0328

<http://www.thecalth.com> info@thecalth.com

MT Promedt Consulting GmbH

Altenhofstrasse 80, 66386 St. Ingbert, Germany
Tel. +49-6894-581020 Fax. +49-6894-581021

SYMBOLS USED

IVD	For in Vitro Diagnostics Use	Temperature limitation "Store at 1-30 °C"
Do not reuse	Used by	
Consult instructions for use	Part number	
LOT	Lot Number	Manufacturer
Sufficient for	Authorized representative in the European Community	
Caution	CE marking	